

KSS DEMOFILTER Test and Demonstration Unit

Compact Ultrafiltration and Microfiltration Test Unit

KSS DEMOFILTER

A portable filtration unit that is ideal for demonstrations since it produces permeate samples immediately. This unit can be used for small scale production applications or the concentration of small process volumes.

FEATURES

- Easy to assemble and operate
- Produces permeate instantly
- Portable
- Light, weighs just 13 lbs
- Compact design, 10" x 10" x 18"
- Uses standard grounded 110V outlet
- Uses KSS 1" ROMICON[®] hollow fiber cartridges and SUPER-COR[®] test cartridges

KSS DEMOFILTER INCLUDES

- DEMOFILTER unit
- Six bag fillers (100 micron)
- Nylon carrying case

TEST CARTRIDGES and REPLACEMENT PARTS

1" ROMICON HOLLOW FIBER CARTRIDGES


PM5	43 mil	KPN 0720041
PM10	43 mil	KPN 0720034
PM10	60 mil	KPN 0720053
PM30	43 mil	KPN 0720039
PM50	43 mil	KPN 0720043
PM50	60 mil	KPN 0720055
PM50	75 mil	KPN 0720040
PM100	20 mil	KPN 0720047
PM100	43 mil	KPN 0720036
PM500	43 mil	KPN 0720045
PM500	75 mil	KPN 0720026
PM500	106 mil	KPN 0720199
WINEFILTER™ 1018	54 mil	KPN 0720061
XM50	45 mil	KPN 0720052
PURON [®] MP 1"		KPN 1015511

1" CARTRIDGES WITH SINGLE ½" SUPER-COR TUBES

HFM-180	KPN 0720027
HFM-513	KPN 0720028

OTHER PARTS

DEMOFILTER Unit	KPN 0995050
Nylon carrying case	KPN 0618003
Bag filter	KPN 0131445



The information contained in this publication is believed to be accurate and reliable, but is not to be construed as implying any warranty or guarantee of performance. We assume no responsibility, obligation or liability for results obtained or damages incurred through the application of the information contained herein. Refer to Standard Terms and Conditions of Sale and Performance Warranty documentation for additional information.

Koch Separation Solutions, Inc. 850 Main Street, Wilmington, MA 01887
Main: +1-978-694-7000 • Fax: +1-978-657-5208 • Toll Free: +1-888-677-5624

For complete contact information and listing of our global locations, visit www.kochseparation.com

©2020 Koch Separation Solutions, Inc. All rights reserved worldwide. For related trademark information, visit www.kochseparation.com/legal.